Zen Oval ZOCP6636 66" x 36" x 24"

PRODUCT NUMBER: ZOCP6636

TECHNICAL SPECIFICATION	S	
Drain	End drain	
Pump Location(s): Select 3	Std. = Standard / Alt LF/RF Optional	
Bath Weight (empty)	165 lbs	
Gallons to Operate	66-84	
ELECTRICAL SPECIFICATIONS		
Motor/Pump	(3) 115v 20AMP GFCI Circuits	
Inline Stainless Steel Heater	(1) 115v 20AMP GFCI Circuit	

STANDARD FEATURES & ACCESSORIES for COMBINATION PLATINUM SERIES BATHS

JETS

33 Air & Water Massage Jets Back Massage System(s) (7 Jets) Neck Massage Pillow (2 Jets) $(\oplus \oplus)$ Sole Soother Jet(s) 🛞 Accumassage Back & Neck Jets Ultra Flow Jets -Classic Air Jets O EZ Clean Air Jets (optional) X

PUMPS

Two 2.25 Hydropower 3-speed Pumps (End Drain Tubs) 1.25 Air Power Variable Speed, Pulse & Wave heated Air Pump Quiet Pad Anti-Vibration Mounted Pump

HEATERS

Stainless Steel Heater Warm Back Rest Tubing Zone

CONTROLS

Air Switch On/Off for Chromatherapy Lighting Air Switch On/Off with digital Temperature Readout Hydro Pumps Electronic Variable Speed Air Massage Control with Automatic Purge & Dry Cycle Control Air Pump 3 Speed Massage: High, Medium or Low

PURITY

HydroClean Water Purity System HydroCleanse

OPTIONAL FEATURES & ACCESSORIES (Please type in the box below)

COMBINATION PLATINUM SERIES



THERS: 1	
STALLATION	
DECK (D)	
UNDER MOUNT (U)	
ALCOVE (A)	
FREESTANDING (F)	
CLAWFOOT (L)	
PEDESTAL (P)	
	DECK (D) UNDER MOUNT (U) ALCOVE (A) FREESTANDING (F) CLAWFOOT (L)



STANDARD ACCESSORIES

Grab Bar(s) - 2 on End Drain baths - White, Biscuit, Bone,
Brushed Nickel, Polished Chrome
Chromatherapy Lighting
Deluxe Cable Drain - White, Biscuit, Bone, Brushed Nickel,
Oil Rubbed Bronze, Polished Chrome
Deluxe Cushion Pillow on 2 person baths
Neck Massage Pillow
EZ Level Base
FDA Recipe Non-Toxic Food & Beverage Grade Plumbing (primarily)
30 Point Quality Assurance Factory Test

Jet & Control placement may vary


REQUIREMENTS

Pump Locations: Select 3 Electrical: Four 115v 20AMP GFCI Circuits Cable Drain and Grab Bars: Select Color / Finish Tub should not be supported by the lip in any way. Level tub & screw or glue EZ Level adjustable legs to sub floor. CUT OUT DIMENSIONS: OVALS: A template will be provided

RECTANGLES: Subtract 1" to 1.5" from length & width for hole cut out size

EZ CLEAN CAST ACRYLIC

- A K HTT® Leg Therapy Air Massage System includes 10 Air Jets on the floor. Optional on Combo Plat and Combo Gold Series
- **B** X 6 Jet HTT® Back Massage System: N/A on Combo Series-Standard on Plat series
- C SW Warm Back Tubing Zone



Some standard & most alternate pump locations may stickout past perimeter of the bath 1"-3". NOTE: Graphics may not be to scale. Specifications and features are subject to change without notice and are ± .25". All final measurements should be taken from tub on site. Access 18"x 18" or larger is required for each pump. Consult your contractor for installation requirements.